

## Product datasheet for **TB415453**

### SCN3A CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of SCN3A (NM_006922), transcript variant 1, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC215453
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<a href="#">NM_006922</a> , <a href="#">NP_008853</a>
Synonyms:	DEE62; EIEE62; FFEVF4; NAC3; Nav1.3
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_008853</a>
Locus ID:	6328
Cytogenetics:	2q24.3
Protein Families:	Druggable Genome, Ion Channels: Sodium, Transmembrane



[View online »](#)