

Product datasheet for **TB415545**

LAP2 (TMPO) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of TMPO (NM_003276), transcript variant 1, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC215545
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	NM_003276 , NP_003267
Synonyms:	CMD1T; LAP2; LEMD4; PRO0868; TP
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_003267
Locus ID:	7112
Cytogenetics:	12q23.1
Protein Families:	Stem cell - Pluripotency, Transmembrane



[View online »](#)