

## Product datasheet for **TB415552**

### ASH2L CytoSection

#### Product data:

|                                       |   |
|---------------------------------------|---|
| Product Type:                         | CytoSections  |
| Description:                          | Transient overexpression of ASH2L (NM_004674), transcript variant 1, in HEK293T cells, paraffin embedded controls for ICC/IHC staining  |
| Species:                              | Human   |
| Expression Host:                      | HEK293T   |
| Expression cDNA Clone or AA Sequence: | TrueORF Clone RC215552  |
| Tag:                                  | C-MYC/DDK   |
| Detection Antibodies:                 | DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)  |
| Target Detection Antibodies:          | ASH2L Mouse Monoclonal Antibody [Clone ID: OTI2H3] (TA810988)   |
| ACCN:                                 | <a href="#">NM_004674</a> , <a href="#">NP_004665</a>   |
| Synonyms:                             | ASH2; ASH2L1; ASH2L2; Bre2  |
| Storage:                              | Room Temperature, or 2-8°C for long term storage  |
| Stability:                            | Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.   |
| Preparation:                          | HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. |
| Note:                                 | This product is for research use only and is not approved for use in humans or in clinical diagnosis.   |
| RefSeq:                               | <a href="#">NP_004665</a>   |
| Locus ID:                             | 9070  |
| Cytogenetics:                         | 8p11.23   |
| Protein Families:                     | Druggable Genome, Transcription Factors   |



[View online »](#)