

Product datasheet for **TB415716**

KCNC1 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of KCNC1 (NM_004976), transcript variant B, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC215716
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	NM_004976 , NP_004967
Synonyms:	EPM7; KV3.1; KV4; NGK2
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_004967
Locus ID:	3746
Cytogenetics:	11p15.1
Protein Families:	Druggable Genome, Ion Channels: Potassium, Transmembrane



[View online »](#)