

## Product datasheet for **TB416035**

### BMP1 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of BMP1 (NM_006128), transcript variant BMP1-2, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC216035
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	BMP1 Mouse Monoclonal Antibody [Clone ID: OT13E9] (TA505813)
ACCN:	<a href="#">NM_006128</a> , <a href="#">NP_006119</a>
Synonyms:	bone morphogenetic protein 1; FLJ44432; PCOLC; PCP; pCP-2; procollagen C-endopeptidase; procollagen C-proteinase; TLD
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_006119</a>
Locus ID:	649
Cytogenetics:	8p21.3
Protein Families:	Druggable Genome, Protease



[View online »](#)