

Product datasheet for **TB416151**

Von Hippel Lindau (VHL) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of VHL (NM_000551), transcript variant 1, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC216151
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	Von Hippel Lindau (VHL) Mouse Monoclonal Antibody [Clone ID: OT11E1] (TA506222)
ACCN:	NM_000551 , NP_000542
Synonyms:	HRCA1; pVHL; RCA1; VHL1
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_000542
Locus ID:	7428
Cytogenetics:	3p25.3
Protein Families:	Druggable Genome, Transcription Factors
Protein Pathways:	Pathways in cancer, Renal cell carcinoma, Ubiquitin mediated proteolysis



[View online »](#)