

## Product datasheet for **TB416155**

### IL32 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of IL32 (NM_001012636), transcript variant 7, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC216155
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<u><a href="#">NM_001012636</a></u> , <u><a href="#">NP_001012654</a></u>
Synonyms:	IL-32alpha; IL-32beta; IL-32delta; IL-32gamma; NK4; TAIF; TAIFa; TAIFb; TAIFc; TAIFd
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u><a href="#">NP_001012654</a></u>
Locus ID:	9235
Cytogenetics:	16p13.3
Protein Families:	Secreted Protein



[View online »](#)