

Product datasheet for **TB416174**

MSMB CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of MSMB (NM_138634), transcript variant PSP57, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC216174
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	MSMB Mouse Monoclonal Antibody [Clone ID: OTI6C7] (TA501072)
ACCN:	NM_138634 , NP_619540
Synonyms:	HPC13; IGBF; MSP; MSPB; PN44; PRPS; PSP; PSP-94; PSP57; PSP94
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_619540
Locus ID:	4477
Cytogenetics:	10q11.22
Protein Families:	Secreted Protein, Transmembrane



[View online »](#)