

Product datasheet for **TB416278**

CYP26A1 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of CYP26A1 (NM_000783), transcript variant 1, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC216278
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	NM_000783 , NP_000774
Synonyms:	CP26; CYP26; P450RAI; P450RAI1
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_000774
Locus ID:	1592
Cytogenetics:	10q23.33
Protein Families:	Druggable Genome, P450, Transmembrane
Protein Pathways:	Retinol metabolism



[View online »](#)