

Product datasheet for **TB416320**

SIRP alpha (SIRPA) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of SIRPA (NM_001040023), transcript variant 2, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC216320
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	SIRP alpha (SIRPA) Mouse Monoclonal Antibody [Clone ID: OTI7B3] (TA807751)
ACCN:	NM_001040023 , NP_001035112
Synonyms:	BIT; CD172A; MFR; MYD-1; P84; PTPNS1; SHPS1; SIRP
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001035112
Locus ID:	140885
Cytogenetics:	20p13
Protein Families:	Druggable Genome, Phosphatase, Transmembrane



[View online »](#)