

## Product datasheet for **TB416326**

### DGKZ CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of DGKZ (NM_201533), transcript variant 3, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC216326
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	DGKZ Mouse Monoclonal Antibody [Clone ID: OT11A6] (TA809946)
ACCN:	<a href="#">NM_201533</a> , <a href="#">NP_963291</a>
Synonyms:	DAGK5; DAGK6; DGK-ZETA; hDGKzeta
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_963291</a>
Locus ID:	8525
Cytogenetics:	11p11.2
Protein Families:	Druggable Genome
Protein Pathways:	Glycerolipid metabolism, Glycerophospholipid metabolism, Metabolic pathways, Phosphatidylinositol signaling system



[View online »](#)