

Product datasheet for **TB416665**

MASP2 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of MASP2 (NM_006610), transcript variant 1, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC216665
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	MASP2 Mouse Monoclonal Antibody [Clone ID: OT12F10] (TA812533)
ACCN:	NM_006610 , NP_006601
Synonyms:	MAP-2; MAP19; MASP-2; MASP1P1; sMAP
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_006601
Locus ID:	10747
Cytogenetics:	1p36.22
Protein Families:	Druggable Genome, Protease, Secreted Protein
Protein Pathways:	Complement and coagulation cascades



[View online »](#)