

Product datasheet for **TB417107**

POFUT2 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of POFUT2 (NM_015227), transcript variant 1, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC217107
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	POFUT2 Mouse Monoclonal Antibody [Clone ID: OT13A12] (TA504972)
ACCN:	NM_015227 , NP_056042
Synonyms:	C21orf80; FUT13
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_056042
Locus ID:	23275
Cytogenetics:	21q22.3



[View online »](#)