

Product datasheet for **TB417302**

HMGA1 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of HMGA1 (NM_145904), transcript variant 6, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC217302
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<u>NM_145904</u> , <u>NP_665911</u>
Synonyms:	AL023995; high mobility group AT-hook 1; high mobility group protein I; HMG-I(Y); HMG-R, HMGIY, HMGA1A, MGC4242, MGC4854, MGC12816; Hmga1a; Hmga1b; Hmgi; HMGI(Y); Hmgiy; HMGY; MGC102580
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_665911</u>
Locus ID:	3159
Cytogenetics:	6p21.31
Protein Families:	Druggable Genome, Stem cell - Pluripotency, Stem cell relevant signaling - JAK/STAT signaling pathway, Transcription Factors



[View online »](#)