

## Product datasheet for **TB417439**

### **NPFF1 Receptor (NPFFR1) CytoSection**

#### **Product data:**

<b>Product Type:</b>	CytoSections
<b>Description:</b>	Transient overexpression of NPFFR1 (NM_022146) in HEK293T cells paraffin embedded controls for ICC/IHC staining
<b>Species:</b>	Human
<b>Expression Host:</b>	HEK293T
<b>Expression cDNA Clone or AA Sequence:</b>	TrueORF Clone RC217439
<b>Tag:</b>	C-MYC/DDK
<b>Detection Antibodies:</b>	Clone OT14C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
<b>ACCN:</b>	<a href="#">NM_022146</a> , <a href="#">NP_071429</a>
<b>Synonyms:</b>	GPR147; NPFF1; NPFF1R1; OT7T022
<b>Storage:</b>	Room Temperature, or 2-8°C for long term storage
<b>Stability:</b>	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
<b>Preparation:</b>	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
<b>Note:</b>	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
<b>RefSeq:</b>	<a href="#">NP_071429</a>
<b>Locus ID:</b>	64106
<b>Cytogenetics:</b>	10q22.1
<b>Protein Families:</b>	Druggable Genome, Transmembrane
<b>Protein Pathways:</b>	Neuroactive ligand-receptor interaction



[View online »](#)