

Product datasheet for **TB417960**

RBFOX1 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of RBFOX1 (NM_018723), transcript variant 4, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC217960
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	RBFOX1 Mouse Monoclonal Antibody [Clone ID: OTI4F9] (TA803488)
ACCN:	NM_018723 , NP_061193
Synonyms:	2BP1; A2BP1; FOX-1; FOX1; HRNBP1
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_061193
Locus ID:	54715
Cytogenetics:	16p13.3



[View online »](#)