

## Product datasheet for **TB418018**

### SCN2A CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of SCN2A (NM_001040143), transcript variant 3, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC218018
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<a href="#">NM_001040143</a> , <a href="#">NP_001035233</a>
Synonyms:	BFIC3; BFIS3; BFNIS; DEE11; EA9; EIEE11; HBA; HBSCI; HBSCII; Na(v)1.2; NAC2; Nav1.2; SCN2A1; SCN2A2
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_001035233</a>
Locus ID:	6326
Cytogenetics:	2q24.3
Protein Families:	Druggable Genome, Ion Channels: Sodium, Transmembrane



[View online »](#)