

## Product datasheet for **TB418549**

### **IQGAP1 CytoSection**

#### **Product data:**

<b>Product Type:</b>	CytoSections
<b>Description:</b>	Transient overexpression of IQGAP1 (NM_003870) in HEK293T cells paraffin embedded controls for ICC/IHC staining
<b>Species:</b>	Human
<b>Expression Host:</b>	HEK293T
<b>Expression cDNA Clone or AA Sequence:</b>	TrueORF Clone RC218549
<b>Tag:</b>	C-MYC/DDK
<b>Detection Antibodies:</b>	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
<b>Target Detection Antibodies:</b>	IQGAP1 Mouse Monoclonal Antibody [Clone ID: OTI6E1] (TA809921)
<b>ACCN:</b>	<a href="#">NM_003870</a> , <a href="#">NP_003861</a>
<b>Synonyms:</b>	HUMORFA01; p195; SAR1
<b>Storage:</b>	Room Temperature, or 2-8°C for long term storage
<b>Stability:</b>	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
<b>Preparation:</b>	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
<b>Note:</b>	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
<b>RefSeq:</b>	<a href="#">NP_003861</a>
<b>Locus ID:</b>	8826
<b>Cytogenetics:</b>	15q26.1
<b>Protein Families:</b>	Druggable Genome
<b>Protein Pathways:</b>	Adherens junction, Regulation of actin cytoskeleton



[View online »](#)