

## Product datasheet for **TB418606**

### Asialoglycoprotein Receptor 2 (ASGR2) CytoSection

#### Product data:

|                                       |   |
|---------------------------------------|---|
| Product Type:                         | CytoSections  |
| Description:                          | Transient overexpression of ASGR2 (NM_080912), transcript variant H2', in HEK293T cells, paraffin embedded controls for ICC/IHC staining  |
| Species:                              | Human   |
| Expression Host:                      | HEK293T   |
| Expression cDNA Clone or AA Sequence: | TrueORF Clone RC218606  |
| Tag:                                  | C-MYC/DDK   |
| Detection Antibodies:                 | DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)  |
| Target Detection Antibodies:          | Asialoglycoprotein Receptor 2 (ASGR2) Mouse Monoclonal Antibody [Clone ID: OTI2A12] (TA505099)  |
| ACCN:                                 | <a href="#">NM_080912</a> , <a href="#">NP_550434</a>   |
| Synonyms:                             | ASGP-R2; ASGPR2; CLEC4H2; HBXBP; HL-2   |
| Storage:                              | Room Temperature, or 2-8°C for long term storage  |
| Stability:                            | Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.   |
| Preparation:                          | HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. |
| Note:                                 | This product is for research use only and is not approved for use in humans or in clinical diagnosis.   |
| RefSeq:                               | <a href="#">NP_550434</a>   |
| Locus ID:                             | 433   |
| Cytogenetics:                         | 17p13.1   |
| Protein Families:                     | Druggable Genome, Transmembrane   |



[View online »](#)