

## Product datasheet for **TB419079**

### ALAS2 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of ALAS2 (NM_000032), transcript variant 1, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC219079
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	ALAS2 Mouse Monoclonal Antibody [Clone ID: OTI2D3] (TA809959)
ACCN:	<u><a href="#">NM_000032</a></u> , <u><a href="#">NP_000023</a></u>
Synonyms:	ALAS-E; ALASE; ANH1; ASB; SIDBA1; XLDPP; XLEPP; XLSA
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u><a href="#">NP_000023</a></u>
Locus ID:	212
Cytogenetics:	Xp11.21
Protein Families:	Druggable Genome
Protein Pathways:	Glycine, serine and threonine metabolism, Metabolic pathways, Porphyrin and chlorophyll metabolism



[View online »](#)