

Product datasheet for **TB419240**

hnRNP L (HNRNPL) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of HNRNPL (NM_001533), transcript variant 1, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC219240
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	hnRNP L (HNRNPL) Mouse Monoclonal Antibody [Clone ID: OT19G7] (TA805445)
ACCN:	<u>NM_001533</u> , <u>NP_001524</u>
Synonyms:	hnRNP-L; HNRPL; P/OKcl.14
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_001524</u>
Locus ID:	3191
Cytogenetics:	19q13.2



[View online »](#)