

Product datasheet for **TB419331**

CSPS (SULT1A3) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of SULT1A3 (NM_001017387), transcript variant 3, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC219331
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	CSPS (SULT1A3) Mouse Monoclonal Antibody [Clone ID: OTI3G5] (TA811313)
ACCN:	NM_001017387 , NP_001017387
Synonyms:	STM, HAST, HAST3, M-PST, ST1A5, TL-PST, SULT1A4, MGC117469
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001017387
Locus ID:	6818
Cytogenetics:	16p11.2
Protein Pathways:	Sulfur metabolism



[View online »](#)