

Product datasheet for **TB419641**

SCDGFB (PDGFD) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of PDGFD (NM_025208), transcript variant 1, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC219641
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	SCDGFB (PDGFD) Mouse Monoclonal Antibody [Clone ID: OTI1C1] (TA503700)
ACCN:	NM_025208 , NP_079484
Synonyms:	IEGF; MSTP036; SCDGF-B; SCDGFB
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_079484
Locus ID:	80310
Cytogenetics:	11q22.3
Protein Families:	Druggable Genome, Secreted Protein
Protein Pathways:	Focal adhesion, Gap junction, Melanoma, Prostate cancer, Regulation of actin cytoskeleton



[View online »](#)