

Product datasheet for **TB419910**

EIF4G1 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of EIF4G1 (NM_198242), transcript variant 4, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC219910
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	NM_198242 , NP_937885
Synonyms:	EIF-4G1; EIF4F; EIF4G; EIF4GI; P220; PARK18
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_937885
Locus ID:	1981
Cytogenetics:	3q27.1
Protein Pathways:	Viral myocarditis



[View online »](#)