

## Product datasheet for **TB420000**

### PCSK9 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of PCSK9 (NM_174936) in HEK293T cells paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC220000
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
Target Detection Antibodies:	PCSK9 Mouse Monoclonal Antibody [Clone ID: OTI4G9] (TA811839)
ACCN:	<a href="#">NM_174936</a> , <a href="#">NP_777596</a>
Synonyms:	FH3; FHCL3; HCHOLA3; LDLCQ1; NARC-1; NARC1; PC9
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_777596</a>
Locus ID:	255738
Cytogenetics:	1p32.3
Protein Families:	Secreted Protein



[View online »](#)