

Product datasheet for **TB420459**

Natriuretic Peptide Receptor B (NPR2) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of NPR2 (NM_003995) in HEK293T cells paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC220459
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
ACCN:	NM_003995 , NP_003986
Synonyms:	AMDM; ANPb; ANPRB; ECDM; GC-B; GCB; GUC2B; GUCY2B; NPRB; NPRBi; SNSK
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_003986
Locus ID:	4882
Cytogenetics:	9p13.3
Protein Families:	Druggable Genome, Protein Kinase
Protein Pathways:	Purine metabolism, Vascular smooth muscle contraction



[View online »](#)