

## Product datasheet for **TB420789**

### CACNG5 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of CACNG5 (NM_014404), transcript variant 2, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC220789
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<u>NM_014404</u> , <u>NP_055219</u>
Synonyms:	calcium channel, voltage-dependent, gamma subunit 5; MGC126656; MGC126682; neuronal voltage-gated calcium channel gamma-5 subunit; voltage-dependent calcium channel gamma-5 subunit
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_055219</u>
Locus ID:	27091
Cytogenetics:	17q24.2
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	Arrhythmogenic right ventricular cardiomyopathy (ARVC), Cardiac muscle contraction, Dilated cardiomyopathy, Hypertrophic cardiomyopathy (HCM), MAPK signaling pathway



[View online »](#)