

Product datasheet for **TB421006**

MTF2 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of MTF2 (NM_007358), transcript variant 1, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC221006
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	NM_007358 , NP_031384
Synonyms:	dj976O13.2; M96; PCL2; TDRD19A
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_031384
Locus ID:	22823
Cytogenetics:	1p22.1
Protein Families:	Druggable Genome, Transcription Factors



[View online »](#)