

Product datasheet for **TB421341**

AGAP1 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of AGAP1 (NM_001037131), transcript variant 1, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC221341
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	AGAP1 Mouse Monoclonal Antibody [Clone ID: OT14E8] (TA504204)
ACCN:	NM_001037131 , NP_001032208
Synonyms:	AGAP-1; CENTG2; cnt-g2; GGAP1
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001032208
Locus ID:	116987
Cytogenetics:	2q37.2
Protein Pathways:	Endocytosis



[View online »](#)