

Product datasheet for **TB421446**

alpha Synuclein (SNCA) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of SNCA (NM_007308), transcript variant 4, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC221446
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	alpha Synuclein (SNCA) Mouse Monoclonal Antibody [Clone ID: OT12F6] (TA506533)
ACCN:	<u>NM_007308</u> , <u>NP_009292</u>
Synonyms:	NACP; PARK1; PARK4; PD1
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_009292</u>
Locus ID:	6622
Cytogenetics:	4q22.1
Protein Families:	Druggable Genome
Protein Pathways:	Alzheimer's disease, Parkinson's disease



[View online »](#)