

## Product datasheet for **TB421594**

### TAF1L CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of TAF1L (NM_153809) in HEK293T cells paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC221594
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
ACCN:	<u><a href="#">NM_153809</a></u> , <u><a href="#">NP_722516</a></u>
Synonyms:	TAF(II)210; TAF2A2
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u><a href="#">NP_722516</a></u>
Locus ID:	138474
Cytogenetics:	9p21.1
Protein Families:	Protein Kinase, Transcription Factors
Protein Pathways:	Basal transcription factors



[View online »](#)