

## Product datasheet for **TB421637**

### MAP2 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of MAP2 (NM_031845), transcript variant 2, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC221637
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	MAP2 Mouse Monoclonal Antibody [Clone ID: OT15D1] (TA813550)
ACCN:	<a href="#">NM_031845</a> , <a href="#">NP_114033</a>
Synonyms:	MAP-2; MAP2A; MAP2B; MAP2C
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_114033</a>
Locus ID:	4133
Cytogenetics:	2q34
Protein Families:	Adult stem cells, Druggable Genome, ES Cell Differentiation/IPS



[View online »](#)