

Product datasheet for **TB421836**

Optineurin (OPTN) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of OPTN (NM_001008213), transcript variant 4, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC221836
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	Optineurin (OPTN) Mouse Monoclonal Antibody [Clone ID: OTI7H2] (TA812132)
ACCN:	NM_001008213 , NP_001008214
Synonyms:	ALS12; FIP2; GLC1E; HIP7; HYPL; NRP; TFIIIA-INTP
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001008214
Locus ID:	10133
Cytogenetics:	10p13
Protein Families:	Druggable Genome



[View online »](#)