

## Product datasheet for **TB421951**

### MLC1 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of MLC1 (NM_139202), transcript variant 2, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC221951
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<u><a href="#">NM_139202</a></u> , <u><a href="#">NP_631941</a></u>
Synonyms:	LVM; MLC; VL
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u><a href="#">NP_631941</a></u>
Locus ID:	23209
Cytogenetics:	22q13.33
Protein Families:	Ion Channels: Other, Transmembrane



[View online »](#)