

## Product datasheet for **TB422041**

### FBXO11 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of FBXO11 (NM_025133), transcript variant 1, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC222041
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	FBXO11 Mouse Monoclonal Antibody [Clone ID: OTI1F6] (TA802323)
ACCN:	<a href="#">NM_025133</a> , <a href="#">NP_079409</a>
Synonyms:	FBX11; IDDFBA; PRMT9; UBR6; UG063H01; VIT1
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_079409</a>
Locus ID:	80204
Cytogenetics:	2p16.3
Protein Families:	Druggable Genome



[View online »](#)