

## Product datasheet for **TB422145**

### KCNH6 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of KCNH6 (NM_173092), transcript variant 2, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC222145
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<u><a href="#">NM_173092</a></u> , <u><a href="#">NP_775115</a></u>
Synonyms:	ERG-2; ERG2; hERG-2; HERG2; Kv11.2
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u><a href="#">NP_775115</a></u>
Locus ID:	81033
Cytogenetics:	17q23.3
Protein Families:	Druggable Genome, Ion Channels: Other, Transmembrane



[View online »](#)