

Product datasheet for **TB422161**

SCN1B CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of SCN1B (NM_199037), transcript variant b, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC222161
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	NM_199037 , NP_950238
Synonyms:	ATFB13; BRGDA5; DEE52; EIEE52; GEFSP1
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_950238
Locus ID:	6324
Cytogenetics:	19q13.11
Protein Families:	Druggable Genome, Ion Channels: Sodium, Transmembrane



[View online »](#)