

Product datasheet for **TB422207**

Desmocollin 2 (DSC2) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of DSC2 (NM_004949), transcript variant Dsc2b, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC222207
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<u>NM_004949</u> , <u>NP_004940</u>
Synonyms:	ARVD11; CDHF2; DG2; DGII/III; DSC3
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_004940</u>
Locus ID:	1824
Cytogenetics:	18q12.1
Protein Families:	Transmembrane
Protein Pathways:	Arrhythmogenic right ventricular cardiomyopathy (ARVC)



[View online »](#)