

## Product datasheet for **TB422718**

### **MAD1 (MAD1L1) CytoSection**

#### **Product data:**

<b>Product Type:</b>	CytoSections
<b>Description:</b>	Transient overexpression of MAD1L1 (NM_001013836), transcript variant 2, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
<b>Species:</b>	Human
<b>Expression Host:</b>	HEK293T
<b>Expression cDNA Clone or AA Sequence:</b>	TrueORF Clone RC222718
<b>Tag:</b>	C-MYC/DDK
<b>Detection Antibodies:</b>	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
<b>ACCN:</b>	<u>NM_001013836</u> , <u>NP_001013858</u>
<b>Synonyms:</b>	MAD1; PIG9; TP53I9; TXBP181
<b>Storage:</b>	Room Temperature, or 2-8°C for long term storage
<b>Stability:</b>	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
<b>Preparation:</b>	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
<b>Note:</b>	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
<b>RefSeq:</b>	<u>NP_001013858</u>
<b>Locus ID:</b>	8379
<b>Cytogenetics:</b>	7p22.3
<b>Protein Families:</b>	Druggable Genome
<b>Protein Pathways:</b>	Cell cycle



[View online »](#)