

Product datasheet for **TB422845**

Antithrombin III (SERPINC1) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of SERPINC1 (NM_000488) in HEK293T cells paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC222845
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
Target Detection Antibodies:	Antithrombin III (SERPINC1) Mouse Monoclonal Antibody [Clone ID: OTI8D5] (TA500455)
ACCN:	NM_000488 , NP_000479
Synonyms:	AT3; AT3D; ATIII; ATIII-R2; ATIII-T1; ATIII-T2; THPH7
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_000479
Locus ID:	462
Cytogenetics:	1q25.1
Protein Families:	Druggable Genome, Secreted Protein
Protein Pathways:	Complement and coagulation cascades



[View online »](#)