

Product datasheet for **TB423032**

MRP4 (ABCC4) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of ABCC4 (NM_005845), transcript variant 1, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC223032
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<u>NM_005845</u> , <u>NP_005836</u>
Synonyms:	MOAT-B; MOATB; MRP4
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_005836</u>
Locus ID:	10257
Cytogenetics:	13q32.1
Protein Families:	Druggable Genome, Ion Channels: Other, Transmembrane
Protein Pathways:	ABC transporters



[View online »](#)