

Product datasheet for **TB423061**

CD10 (MME) CytoSection

Product data:

| | |
|---------------------------------------|---|
| Product Type: | CytoSections |
| Description: | Transient overexpression of MME (NM_007287), transcript variant 1bis, in HEK293T cells, paraffin embedded controls for ICC/IHC staining |
| Species: | Human |
| Expression Host: | HEK293T |
| Expression cDNA Clone or AA Sequence: | TrueORF Clone RC223061 |
| Tag: | C-MYC/DDK |
| Detection Antibodies: | DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569) |
| Target Detection Antibodies: | CD10 (MME) Mouse Monoclonal Antibody [Clone ID: OT13D11] (TA504507) |
| ACCN: | NM_007287 , NP_009218 |
| Synonyms: | CALLA; CD10; CMT2T; NEP; SCA43; SFE |
| Storage: | Room Temperature, or 2-8°C for long term storage |
| Stability: | Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed. |
| Preparation: | HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. |
| Note: | This product is for research use only and is not approved for use in humans or in clinical diagnosis. |
| RefSeq: | NP_009218 |
| Locus ID: | 4311 |
| Cytogenetics: | 3q25.2 |
| Protein Families: | Druggable Genome, Protease, Transmembrane |
| Protein Pathways: | Alzheimer's disease, Hematopoietic cell lineage, Renin-angiotensin system |



[View online »](#)