

## Product datasheet for **TB423116**

### CD10 (MME) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of MME (NM_007289), transcript variant 2b, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC223116
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	CD10 (MME) Mouse Monoclonal Antibody [Clone ID: OT13D11] (TA504507)
ACCN:	<a href="#">NM_007289</a> , <a href="#">NP_009220</a>
Synonyms:	CALLA; CD10; CMT2T; NEP; SCA43; SFE
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_009220</a>
Locus ID:	4311
Cytogenetics:	3q25.2
Protein Families:	Druggable Genome, Protease, Transmembrane
Protein Pathways:	Alzheimer's disease, Hematopoietic cell lineage, Renin-angiotensin system



[View online »](#)