

Product datasheet for **TB423430**

Dystonin (DST) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of DST (NM_015548), transcript variant 1eA, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC223430
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	NM_015548 , NP_056363
Synonyms:	BP240; BPA; BPAG1; CATX-15; CATX15; D6S1101; DMH; DT; EBSB2; HSAN6; MACF2
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_056363
Locus ID:	667
Cytogenetics:	6p12.1



[View online »](#)