

Product datasheet for TB423522

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

UBE2J2 CytoSection

Product data:

Product Type: CytoSections

Description: Transient overexpression of UBE2J2 (NM_058167), transcript variant 2, in HEK293T cells,

paraffin embedded controls for ICC/IHC staining

Species: Human
Expression Host: HEK293T

Expression cDNA Clone

or AA Sequence:

TrueORF Clone RC223522

Tag: C-MYC/DDK

Detection Antibodies: DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)

ACCN: <u>NM 058167</u>, <u>NP 477515</u>

Synonyms: NCUBE-2; NCUBE2; PRO2121

Storage: Room Temperature, or 2-8°C for long term storage

Stability: Blocks are guaranteed for a year from the date of receipt if proper storage instructions were

followed.

Preparation: HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells

were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin &

dehydrated before embedding in paraffin.

Note: This product is for research use only and is not approved for use in humans or in clinical

diagnosis.

RefSeq: NP 477515

Locus ID: 118424 **Cytogenetics:** 1p36.33

Protein Families: Transmembrane

Protein Pathways: Parkinson's disease, Ubiquitin mediated proteolysis

