

Product datasheet for **TB423548**

H2AZ2 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of H2AFV (NM_201436), transcript variant 3, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC223548
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	NM_201436 , NP_958844
Synonyms:	H2A.Z-2; H2AFV; H2AV
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_958844
Locus ID:	94239
Cytogenetics:	7p13
Protein Families:	Druggable Genome
Protein Pathways:	Systemic lupus erythematosus



[View online »](#)