

Product datasheet for **TB423673**

Metabotropic Glutamate Receptor 4 (GRM4) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of GRM4 (NM_000841) in HEK293T cells paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC223673
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
Target Detection Antibodies:	Metabotropic Glutamate Receptor 4 (GRM4) Mouse Monoclonal Antibody [Clone ID: OTI6D1] (TA810384)
ACCN:	NM_000841 , NP_000832
Synonyms:	GPRC1D; mGlu4; MGLUR4
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_000832
Locus ID:	2914
Cytogenetics:	6p21.31
Protein Families:	Druggable Genome, GPCR, Transmembrane
Protein Pathways:	Neuroactive ligand-receptor interaction, Taste transduction



[View online »](#)