

Product datasheet for **TB423704**

CTNND1 CytoSection

Product data:

| | |
|---------------------------------------|---|
| Product Type: | CytoSections |
| Description: | Transient overexpression of CTNND1 (NM_001085469), transcript variant 13, in HEK293T cells, paraffin embedded controls for ICC/IHC staining |
| Species: | Human |
| Expression Host: | HEK293T |
| Expression cDNA Clone or AA Sequence: | TrueORF Clone RC223704 |
| Tag: | C-MYC/DDK |
| Detection Antibodies: | DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569) |
| Target Detection Antibodies: | CTNND1 Mouse Monoclonal Antibody [Clone ID: OT12E8] (TA800807) |
| ACCN: | NM_001085469 , NP_001078938 |
| Synonyms: | BCDS2; CAS; CTNND; p120; p120(CAS); p120(CTN); P120CAS; P120CTN |
| Storage: | Room Temperature, or 2-8°C for long term storage |
| Stability: | Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed. |
| Preparation: | HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. |
| Note: | This product is for research use only and is not approved for use in humans or in clinical diagnosis. |
| RefSeq: | NP_001078938 |
| Locus ID: | 1500 |
| Cytogenetics: | 11q12.1 |
| Protein Families: | Druggable Genome |
| Protein Pathways: | Adherens junction, Leukocyte transendothelial migration |



[View online »](#)