

Product datasheet for **TB423907**

VE Cadherin (CDH5) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of CDH5 (NM_001795) in HEK293T cells paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC223907
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
Target Detection Antibodies:	VE Cadherin (CDH5) Mouse Monoclonal Antibody [Clone ID: OTI1F4] (TA804746)
ACCN:	<u>NM_001795</u> , <u>NP_001786</u>
Synonyms:	7B4; CD144
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_001786</u>
Locus ID:	1003
Cytogenetics:	16q21
Protein Families:	Druggable Genome, ES Cell Differentiation/IPs, Transmembrane
Protein Pathways:	Cell adhesion molecules (CAMs), Leukocyte transendothelial migration



[View online »](#)