

Product datasheet for **TB423951**

SCN4B CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of SCN4B (NM_174934), transcript variant 1, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC223951
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	NM_174934 , NP_777594
Synonyms:	ATFB17; LQT10; Navbeta4
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_777594
Locus ID:	6330
Cytogenetics:	11q23.3
Protein Families:	Ion Channels: Sodium, Transmembrane



[View online »](#)